

LIFTER Ascend



JOMOR
HEALTHCARE
Quality - Value - Design

Nimble, compact and extremely easy to manoeuvre

With its unrivalled range of movement, the Ascend can transfer a patient from a low-seated position to fully extended with ease. Its active lifting motion was specially designed to encourage user participation, which ultimately promotes patient independence and well-being.



Supportive knee pad: The knee pad is adjustable in height (by 152mm or 6 inches) to ensure that the patient feels comfortable and fully supported throughout the lifting process. The pad, which has a molded design to provide additional support, also has a leg restraint for those patients needing additional reassurance.

The Ascend has the option of both a standing and transport sling, providing the carer the option of performing various care duties, such as:

- Toileting
- Standing
- Point to point transfers

In addition to its lightweight aluminium construction the Ascend has a powerful safe working load of 170kgs/ 375lbs/27st.



LIFTER Ascend

Specification	Imperial	Metric
Safe Working Load	375 lbs	170 kg
Maximum Overall Length	43.3"	1100mm
Minimum Overall Length	39.3"	1000mm
Maximum Overall Height	67.3"	1710mm
Minimum Overall Height	42.1"	1070mm
Maximum Height to Attachment Point	65.5"	1665mm
Minimum Height to Attachment Point	38.9"	990mm
Turning Radius	51.9"	1320mm
Legs Open - External Width	39.1"	995mm
Legs Open - Internal Width	34.8"	885mm
Legs Closed - External Width	24.8"	630mm
Legs Closed - Internal Width	20"	510mm
Overall Height of Legs	4.7"	120mm
Ground Clearance	1.4"	35mm
Front Twin Castors	4"	100mm
Rear Braked Castors	4"	100mm

Weights		
Mast, Base & Boom Assembly	97 lbs	44 kg
Power Pack	6.6 lbs	3 kg
Total	103.6 lbs	47 kg
Base Assembly	42.9 lbs	19.5 kg
Mast & Boom Assembly	45.1 lbs	20.5 kg

Specifications may vary due to manufacturing tolerances

Removable foot tray

Removing the foot tray transforms the Ascend into a device that allows a patient to fully stand. This is made possible by the Ascend's extensive lifting range. Once standing, the lift can be removed and the patient can be encouraged to walk independently.



For more information contact:

Jomor Healthcare P/L 308 Lower Dandenong Road
 Mordialloc, VIC 3195
ph: 1300 651 235 **fax:** (03) 9588 2255
email: sales@jomor.com.au **web:** www.jomor.com.au

JOMOR
 HEALTHCARE

Quality - Value - Design

